

T4 DNA Ligase



3.2018



blirt

T4 DNA Ligase

T4 DNA Ligase is an ATP-dependent recombinant enzyme isolated from *Escherichia coli* strain used to clone DNA. **T4 DNA Ligase** catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA. It will join both blunt-ended and cohesive-ended restriction fragments of DNA, as well as repair single-stranded nicks in duplex DNA, RNA or DNA/RNA hybrids.



Applications

- Cloning PCR products
- Cloning restriction fragments
- Joining double-stranded oligonucleotide linkers or adaptors to DNA
- Site-directed mutagenesis
- Nick repair in duplex DNA, RNA or DNA/RNA hybrids
- Self-circularization of linear DNA
- LM PCR methods (Ligation Mediated PCR), e.g. amplified fragment length polymorphism (AFLP)

T4 DNA Ligase

Protocol

1. Add the reaction reagents listed below to a sterile nuclease-free Eppendorf tube placed **on ice in a freezing rack**. The reaction agents should be added in the following order:

Component	Volume
Vector DNA	x μ l (20 – 50 ng)
Insert DNA	y μ l (3 – 10 molar excess over vector)*
10x T4 Ligation Buffer	2 μ l
ATP Solution (25 mM)	0.4 – 0.8 μ l
T4 DNA Ligase (5 U/ μ l)	1 μ l
Nuclease-free water	up to 20 μ l

* A lower ratio will result in a less efficient ligation; a higher ratio will incite multiple insertions. This data is intended for use as a guide only; conditions will vary from reaction to reaction and may need optimization.

2. Mix gently and spin briefly.
3. For cohesive (sticky) ends, incubate at 20-25°C for 0.5-1 hour. For blunt ends incubate at 20-25°C for 1-2 hours.
4. Cool the samples on ice and transform 1-5 μ l of the reaction mixture into 50 μ l *E. coli* competent cells.



Additional information

- The 10x T4 Ligation Buffer and ATP Solution should be thawed and resuspended at room temperature.
- For blunt-end ligations, use higher quantities of both the vector and the insert DNA.
- PEG greatly increases the ligation efficiency of blunt-end DNA ligation. The recommended concentration of PEG in the ligation reaction mixture is 5% (w/v) (high concentration of PEG significantly reduces the transformation efficiency of electrocompetent cells).
- For sticky (cohesive)-end ligations, we recommend heating both the vector and the insert DNA prior to ligation.
- The electrotransformation efficiency may be improved by heat inactivation of the **T4 DNA Ligase** and purification of the DNA by means of a spin column purification method (**EXTRACTME DNA CLEAN-UP** kit or **EXTRACTME DNA GEL-OUT** kit).
- We recommend using a 3–10 molar excess of insert DNA over vector DNA. To calculate the optimal amounts of insert DNA in a ligation reaction, use the following equation:

$$\text{ng of insert} = \frac{\text{ng of vector} \times \text{kb size of insert} \times \text{molar ratio insert:vector}}{\text{kb size of vector}}$$

Example: If using 20 ng of a vector plasmid (4 kb), for a 5:1 molar ratio of insert:vector, you will require the following quantity of 1 kb insert:

$$\frac{20 \times 1 \times 5}{4} = 25 \text{ ng}$$

- The enzyme is inhibited by >200 mM NaCl or KCl concentrations.
- Inactivate enzyme at 65°C for 10 minutes or at 70°C for 5 minutes.

Storage buffer

10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 50 mM KCl, 1 mM DTT, 50% (v/v) glycerol

10x T4 Ligation Buffer

660 mM Tris-HCl (pH 7.5), 100 mM DTT, 100 mM MgCl₂

The 10x T4 Ligation Buffer does not contain ATP, which must be added separately. For most applications, ATP should be added to the reaction to a final concentration of 0.5–1.0 mM. 25 mM ATP Solution is included.

50% PEG Solution

50% (w/v) polyethylene glycol 4000.



Quality control

The product was tested in ligation of *Hind*III-cut lambda DNA with a different amounts of enzyme. The product is free of unspecific DNA nucleases.

Unit definition

One (Weiss) unit of **T4 DNA Ligase** catalyzes the conversion of 1 nmol of ^{32}P from pyrophosphate into Norit-adsorbable material in 30 minutes at 37°C.

One Weiss unit is equivalent an approximately 200 cohesive end units.

T4 DNA Ligase

Components	EN11-050 500 U	EN11-250 2500 U
T4 DNA Ligase 5 U/ μ l	100 μ l	5x 100 μ l
10x T4 Ligation Buffer	200 μ l	5x 200 μ l
ATP Solution (25 mM)	80 μ l	5x 80 μ l
50% PEG Solution	200 μ l	5x 200 μ l

Storage & shipping

Storage conditions

All components should be stored at -20°C .

Shipping conditions

Shipping on dry or blue ice.

 For research use only